





## Metrological characterization of standards – plastic films, which are used for calibration of nano volume spectrophotometers

## Veljko Zarubica <sup>1</sup>, Edvard Kubičela<sup>1</sup>, Boban Zarkov<sup>2</sup> and Milesa Srećković<sup>3</sup>

- <sup>1</sup>Laboratory for calibration and validation, Analysis d.o.o, Japanska 4, Belgrade, Serbia
- <sup>2</sup> Bureau of measures and precious metals, Mike Alasa 14, Belgrade, Serbia
- <sup>3</sup>Faculty of Electrical Engineering, Belgrade, Serbia

**Abstract.** UV-VIS spectrophotometers, including nano drop (volume) spectrophotometers, provide an answer to all challenges in scientific and industrial laboratories, whether for pharmaceutical, chemical or service laboratories, all the way to more complex applications such as quantification of DNA, RNA, and protein per second, using only 1-2 µL of the sample. Certified reference materials are used to check the measurement capabilities of nanodrop spectrophotometers and other analytical instruments. Since these are solutions, which by themselves are not practical in terms of handling, stability, availability, we were thinking of using materials for their calibration, which would have more suitable properties than the solutions. We decided to use plastic films, which are thin enough not to be a problem during manipulation and measurement, which are stable in terms of determining the regular transmittance, which are optically neutral in the wider spectrum of optical radiation, and which are temperature stable and do not change their metrological characteristics during prolonged exposure to light. The paper evaluates all these influential characteristics and presents the budget of measurement uncertainty of plastic films in the UV part of the optical radiation spectrum at wavelengths of 260 nm and 280 nm (Fig.6). The term "metrologically characterized materials" represent equivalent of the materials, with precisely define regular transmittance (e.g., absorption), as well as estimated influential parameters which may appear, thus contributing to the uncertainty of the measurement.

Calibration of nano volume spectrophotometers (Fig.3) is usually performed with some of the standard reference solutions. It is foreseen by the manufacturer of these instruments, that the mentioned solutions are applied for one-time use. After using these solutions, any excess solution is thrown away, which is basically the negative side of these SRMs. Of course, the other negative side of this method of calibration is the high price of the SRM. Our wish was to try to find a solution, namely SRM, such that the price of SRM is acceptable and that there is a possibility of using these SRM for a longer period. In this way, we would be able to calibrate Nano volume spectrophotometers with a completely different type of standard reference materials - foils. Our goal was to characterize these foils metrologically and present them as CRMs, which will enable them to be used for a long period of time, without losing, or in the best case, degrading, their metrological properties, especially in terms of long-term stability, accuracy of reproducing the standard values etc. Before the characterization itself, the mentioned standard reference materials-foils were defined in terms of their composition analysis. We used an IR spectrophotometer with Fourier transform (FT-IR) (Fig. 1) to analyse the composition of the mentioned foils in detail. The results of the analysis are presented in Fig. 2.


The influencing parameters that occur during the measurement of the spectral coefficient of transmittance, i.e. absorption, of optically neutral filters are as follows:

- homogenity, exposure, temperature and atmospheric influences, light exposure stability and the influence of dispersion.
- The influencing parameters, such as, homogenity, light exposure stability and the influence of dispersion were determined by the spectrophotometric method using the spectrophotometric system (Fig. 4), while the temperature influence parameter was determined by the method of the integration sphere (Fig. 5)

## Fig.2 Spectral characteristic filter foil obtained FT-IR FOL 01 ANALYSIS Fig.1 IR spectrophotometer with Fourier transform (FT-IR) Fiber Number: A0305 Manufacturer's Code: 375d 0.5 Trade Name: MIRAFI Supplier: MIRAFI INC. Fiber Number: C0413 Manufacturer's Code: 700X Wavenumbers (cm-1) Spectrum: FOL 01 Region: 3495.26-650.06 Number of sample scans: 32 Correlation Number of background scans: 32 Resolution: 4.000

Synthetic Fibers by Microscope

Fig.3 Nano drop spectrophotometer



Index Match Compound name

143 92.15 OLEFIN

374 91.73 OLEFIN

Sample gain: 1.0 Optical velocity: 0.4747

Aperture: 80.00





Fig. 5 Method of determination homogenity filter foils with method wich used the integration sphere

and cleaning withou



Drive with Sliding Coupling

Table 1. Determined of stability during a period of one year

Measuring of the spectral coefficient of transmittance at a determined stability	FOL1		
during a period of one year	230 nm	260 nm	280 nm
240104	0,2275	0,6425	0,6323
240311	0,2244	0,6383	0,6275
240521	0,2204	0,6371	0,6226
240926	0,2277	0,6321	0,6194
241129	0,2233	0,6267	0,6137
Average τ	0,2247	0,6353	0,6231
u component of measurement uncertainty of stability	0,0070	0,0049	0,0059
u ave	0,0104		

Table 2. Measurement uncertainty, resulting from the inhomogenity of the three

spectrophotometric optical neutral folies				
Measuring of the spectral coefficient of transmittance at a determined wavelength before light exposition	FOL1	FOL2	FOL3	
1 central position	0,2268	0,6425	0,6323	
2 left positions	0,2270	0,6401	0,6320	
3 right positions	0,2273	0,6383	0,6275	
4 up positions	0,2280	0,6371	0,6226	
5 down positions	0,2272	0,6341	0,6286	
τ	0,2270	0,6403	0,6306	
u component of measurement uncertainty due to inhomogenity of filter surface	0,0006	0,0019	0,0025	
U ave	0,0032			

Table 3. Measurement uncertainty, resulting

from infuence of light radiation for FOL1

_			
Measuring of the spectral coefficient of transmittance on light exposure for FOL1	230 nm	260 nm	230 nm
1	0,2277	0,6321	0,6194
2	0,2233	0,6267	0,6137
3	0,2242	0,6393	0,6275
Average t	0,2255	0,6294	0,6166
Standard deviation	0,0024	0,0063	0,0069
Relative standard deviation	1,05%	1,00%	1,12%
u component of measurement uncertainty due to exposure to light	-0,0004	0,0028	0,0031
uave		0,0019	

Table 4. Measuring on temperature at 7,5 °C

Measuring results of spectral FOL 2 FOL 1 FOL 3 coefficient transmittance on 7,5 oC Mean value of 0,2222 0,6322 0,6161 series at 6 measuring Standard 0,00010 0,00011 0,00007 deviation 0,00044 0,00017 0,00012 RSD % Measurement 0,00010 0,00025 0,00007 uncertainty 0,00028 Uav

Table 4. Measuring on temperature at 20,5 °C

0.251

Measuring results of spectral coefficient transmittance on 20,5 oC	FOL 1	FOL 2	FOL 3
Mean value of series at 6 measuring	0,2255	0,6288	0,6155
Standard deviation	0,00004	0,00007	0,00003
RSD %	0,00016	0,00011	0,00005
Measurement uncertainty	0,00009	0,00006	0,00003
Uav		0,00011	

Table 5. Budget of the measuring uncertainty of the measurement of the spectral coefficient of transparency FOL1

	FOL1			
The source of measurement uncertainty	Coefficient of sensitivity	Contribution of standard uncertainty (Type A)	Contribution of standard uncertainty (Type B)	
Component of measurement uncertainty of				
spectrophotometric system (dispersion of results,				
drift, stray light, resolution, fault radiation, setting				
wavelength)	1		0,00104	
Long stability of the results of measurement of				
spectral coefficient of transmittance after one year	1	0,0104		
Component of measurement uncertainty of		,		
homogeneity of filter structure	1	0,0006		
Component of measurement uncertainty of exposure				
of light radiation	1	0,0004		
Component of measurement uncertainty of				
temperature influence	1	0,0017		
Standard uncertainty by measurement of spectral		0,0106	0,00104	
transmittance coefficient		0,0100	0,00104	
Combined uncertainty of the measurements of the	0,0213			
transmittance coefficient with k=2				

Table 6. Budget of the measuring uncertainty of the measurement of the spectral coefficient of transparency FOL2

spectral coefficient of transparency FUL2				
The source of measurement uncertainty	FOL2			
	Coefficient of sensitivity	Contribution of standard uncertainty (Type A)	Contribution of standard uncertainty (Type B)	
Component of measurement uncertainty of				
spectrophotometric system (dispersion of results, drift, stray light, resolution, fault radiation, setting				
wavelength)	1		0,00104	
Long stability of the results of measurement of spectral coefficient of transmittance after one year	1	0,0110		
Component of measurement uncertainty of homogeneity of filter structure	1	0,0019		
Component of measurement uncertainty of exposure of light radiation	1	0,0028		
Component of measurement uncertainty of temperature influence	1	0,0002		
Standard uncertainty by measurement of spectral transmittance coefficient		0,0115	0,00104	
Combined uncertainty of the measurements of the transmittance coefficient with k=2		0,0231		

Table 7. Budget of the measuring uncertainty of the measurement of the spectral coefficient of transparency FOL3

The source of measurement uncertainty	FOL3		
	Coefficient of sensitivity	Contribution of standard uncertainty (Type A)	Contribution of standard uncertain (Type B)
Component of measurement uncertainty of spectrophotometric system (dispersion of results, drift, stray light, resolution, fault radiation, setting wavelength)	1		0,00104
Long stability of the results of measurement of spectral coefficient of transmittance after one year	1	0,0173	3,00101
Component of measurement uncertainty of homogeneity of filter structure	1	0,0025	
Component of measurement uncertainty of exposure of light radiation	1	0,0031	
Component of measurement uncertainty of temperature influence	1	0,0001	
Standard uncertainty by measurement of spectral transmittance coefficient		0,0177	0,00104
Combined uncertainty of the measurements of the transmittance coefficient with $k=2$		0,0356	